



scenesafe™

Evidence Recovery Systems

Contents

Why SceneSafe?	2
Product Key & Definitions	3
Forensic Evidence Recovery	4
DNA Sampling	5
Sexual Assault	7
Trace DNA	14
Toxicology.....	17
Firearms	19
Drugs	21
PEK (Personal Exploitation Kit)	24
Index	26





Why SceneSafe?

With our technical expertise, combined experience and worldwide reputation, SceneSafe have developed and stock a full range of products to support the continued advances in Forensic Science.

At SceneSafe we pride ourselves on being at the forefront of forensic development taking technology to new levels. Offering our customers a bespoke solution has been the success of our growing forensic range.

We continue to develop our products taking into account the recommendations of the Faculty of Forensic & Legal Medicine (FFLM) and the UK Forensic Regulator. We also adhere to the Specification for consumables used in the collection, preservation and processing of material for forensic analysis (PAS 377). SceneSafe products are independently validated to a forensic DNA Grade.

SceneSafe have developed the range of products to include four key areas of forensics:

- Evidence Recovery Systems
- Crime Scene Equipment
- Forensic Laboratory Equipment
- Training

By offering the above, we are able to develop a complete forensic capability equipping and training to an industry standard. Our expertise and platform in the market allows us to provide a single-source solution in crime scene programmes.

Our high quality standards have been implemented on the advice of forensic experts and are monitored and reviewed continually. SceneSafe is accredited to ISO 9001 Quality and ISO14001 Environmental standards.

Through our philosophy of total quality, continuous improvement and teamwork, we endeavour to surpass your expectations and provide the highest levels of customer care and satisfaction.

Product Key & Definitions



Protect: Protecting the crime scene and valuable evidence during the early stages of investigation is vital to ensure successful judicial outcomes. SceneSafe has selected the products in this range carefully to maximise evidential integrity whilst enhancing the safety of all personnel involved.



Collect: Collecting evidence quickly with the minimum of contamination will provide robust samples for subsequent examination and forensic analysis. The high quality products identified in this range are practical and ideal for use at crime scenes.



Transport: Maintaining the evidential chain of integrity when moving evidence from crime scene to an examination area is key in police investigations. SceneSafe products in this range are designed using appropriate continuity details and use quality assured tamper evident technology.



Identify: Whether at scene or during a detailed examination, getting results rapidly can add substantial value to any investigation. Products in this SceneSafe range provide quick evaluations of the evidence.

Look out for the specific quality assurances and added value benefits against our SceneSafe products

Forensic DNA Grade:

All critical components used for the detection and collection of DNA material have been quality controlled by testing a representative sample and analysed for the presence of DNA. This process is carried out by ISO 17025 accredited forensic laboratories.

DNA Denatured:

Components are treated with Ethylene Oxide (EtO) to denature any possible background DNA to a point where it cannot be amplified. Ethylene oxide is a gas that kills bacteria, mould and fungi.

ISO17025 Compliant:

Components are validated by an ISO17025 accredited analytical forensic service provider.

Meets FFLM guidelines:

Product follows the recommendations of the Forensic Science Committee of the Faculty of Forensic & Legal Medicine of the Royal College of Physicians. The FFLM is set up to develop and maintain the highest possible standards of competence and professional integrity in forensic and legal medicine.

Drug-free:

All critical components have been batch quality tested by sampling a representative quantity for the presence of drugs and drug metabolites. This is undertaken using best technical methods of analysis.

Alcohol-free:

All critical components have been batch quality tested by sampling a representative quantity for the presence of alcohols. This is undertaken using best technical methods of analysis.

Gunshot Residue (GSR) free:

All critical components have been batch quality tested by sampling a representative quantity for the presence of gunshot residues and other particulates. This is undertaken using best technical methods of analysis.

Produced in accordance with PAS 377:

Each kit is produced following the guidelines provided in the specification for consumables used in the collection, preservation and processing of material for forensic analysis guidelines (PAS 377).

Clean Room Manufactured:

Each kit is manufactured in a ISO 14644 certified class 4 and class 5 clean rooms by staff trained in the production of kits used in the processing / collecting of DNA and other forensic samples. Our clean room protocols include full audit trail of components, environmental monitoring, PPE training, clean room cleaning and batch testing. All production staff DNA profiles are recorded and maintained for elimination purposes. Certification can be provided upon request.

FORENSIC EVIDENCE RECOVERY

scenesafe™
evidence recovery systems

Tel: +44 (0)1621 786654

DNA is unique to each individual in its make up. This means that samples taken containing DNA are of great value when investigating a suspected crime.

At SceneSafe we have been developing evidence recovery kits for many years. Over time we have developed an industry standard method for recovering vital DNA material in a secure and safe manner ensuring integrity is maintained throughout.

In this brochure you will find product information relating to the following areas of DNA recovery:

- DNA Sampling
- Sexual Assault
- Trace DNA
- Toxicology
- Firearms
- Drug Testing

SceneSafe have the capability to design and manufacture custom made DNA evidence recovery kits specific to your requirements. We have worked with Police Forces around the world to develop products to suit their exact needs.

Contact us today and find out how we can support you in the fight against crime.
Tel: +44 (0)1621 786654 Email: Sales@scenesafe.co.uk

Evidence Recovery Systems: DNA Sampling

Every Evidence Recovery Kit is made following a strict quality control procedure following PAS 377 guidelines. With the requirement in today's industry demanding forensic material to be free of contamination we do all that is required to ensure the integrity of products is maintained. The components undergo an external post production validation process so that our kits are certified to the industry standard.



PACE DNA Sampling Kit



- Forensic DNA grade
- DNA Denatured
- ISO17025 Compliant
- Clean Room Manufactured

Enables mouth swabs or plucked hairs to be taken from persons arrested and detained in police custody in respect of a recordable offence. The DNA profile obtained from these samples will be added to the National DNA Database

Code K505 / Pack Size 25



Elimination DNA Sampling Kit



- Forensic DNA grade
- DNA Denatured
- ISO17025 Compliant
- Clean Room Manufactured

Enables mouth swabs to be taken from persons who are not suspects, but who wish to volunteer to provide samples to assist in police investigations. Mainly used for Intelligence-Led Screens.

Code K515 / Pack Size 25



Police Elimination DNA Sampling kit



- Forensic DNA grade
- DNA Denatured
- ISO17025 Compliant
- Clean Room Manufactured

Enables mouth swabs to be taken from Police and support staff for inclusion on the UK Police Elimination DNA Database

Code K530 / Pack Size 25

Evidence Recovery Systems: DNA Sampling

The kits below are examples of the DNA sampling kits designed for use by the Scottish Police Authority (SPA) and authorised to be stored on the Scottish DNA Database. SceneSafe offer custom made DNA sampling kits to collect the data to be later analysed and stored on a national or regional DNA database system.



DNA Arrest/Detention Sampling Kit



- Forensic DNA grade
- DNA Denatured
- ISO17025 Compliant
- Clean Room Manufactured

Enables mouth swabs or hair samples to be collected from persons arrested to detained in police custody in respect of a recordable offence. The DNA profile obtained from these samples will be stored on the Scottish DNA Database as well as the UK NDNAD and can only be used for intelligence purposes and are not normally used in a court of law.

Code K501 / Pack Size 25



DNA Intelligence-Led Screen Sampling Kit



- Forensic DNA grade
- DNA Denatured
- ISO17025 Compliant
- Clean Room Manufactured

Enables collection of buccal samples to be taken from persons who are not suspects, but who wish to volunteer to provide samples to assist in police investigations. Mainly used for Intelligence-Led Screens

Code K511 / Pack Size 25



Police Elimination DNA Sampling Kit



- Forensic DNA grade
- DNA Denatured
- ISO17025 Compliant
- Clean Room Manufactured

Enables mouth swabs to be taken from police and support staff for inclusion on the Scottish Police Elimination DNA Database. Profiles are not held on the Scottish DNA Database.

Code K532 / Pack Size 25

Custom Made Options:

SceneSafe offer custom made DNA Sampling kits designed in any language. Each kit can be tailored to the exact specifications to suit the sampling methods and protocols in place.

Validation:

All of our kits are validated, certified and batch controlled to offer our customers a full traceable audit trail.

Specifications include:

- Sequential numbering
- Tamper evident packaging
- Unique barcode
- Custom logo
- Language option
- DNA forensic grade material

Sexual Assault: Early Evidence

The collection of evidence needs to begin from the moment a sexual assault claim is reported. This is essential to prevent loss or deterioration occurring within the passage of time. SceneSafe has developed a range of modular examination kits in line with the FFLM guidelines and PAS 377.

We can offer custom designed sexual assault kits tailored to your exact specifications. Contact a member of our development team for more information.

Our standard kits are available in stock ready for deployment to any location in the world.



Evidence Recovery Systems: Sexual Assault



Early Evidence Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

A first response kit for the collection of samples at the earliest opportunity to capture vital evidence. An early urine sample may increase the chances of detecting drugs and/or alcohol that may have been used to facilitate a sexual assault. A mouth rinse sample may detect vital evidence where oral penetration has been alleged. In addition it avoids unnecessary discomfort for the complainant whilst a full medical examination is awaited.

Code K106 / Pack Size 1



**Early Evidence Kit
with swabs**



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

The addition of swabs to the standard early evidence kit allows for sampling of areas of the skin where body fluid traces may have been deposited e.g. bites. It is recommended that mouth swabs are only taken once a medical examiner has had a chance to assess the mouth of the complainant for damage/injury.

Code K106SW / Pack Size 1



**Early Evidence Kit
(Two urine)**



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

This kit contains two urine toilet paper modules as a second urine sample is recommended where the incident took place in the preceding 24 hours.

Code K108 / Pack Size 1



**Early Evidence Kit
(Two urine) with swabs**



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

As with the K108 this kit features two urine toilet paper modules with the addition of skin swabs, this is a more comprehensive kit which could travel with the complainant to take any necessary early evidence samples in the time leading up to the medical examination.

Code K108SW / Pack Size 1

Evidence Recovery Systems: Sexual Assault



Alcohol/Drug Blood Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

Alcohol/Drug Urine Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

To collect and store a urine sample suitable for subsequent toxicological analysis to detect alcohol, drugs of abuse, medicinal and sedative drugs, including drugs that may be used in a drug-facilitated sexual assault.

Mouth Collection Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

To collect and store samples from the mouth of a complainant suitable for subsequent forensic analysis to detect semen and/or DNA. Contains swabs and return container.

Standard Swab Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

Designed for the recovery of body fluids or lubricant from the complainant or suspect. The components are suitable for samples requiring analysis using the latest DNA profiling techniques. This kit is also available with an

Code K380 / Pack Size 1

Code K381 / Pack Size 1

Code K383 / Pack Size 1

Code K392 / Pack Size 1

Code K392W / Pack Size 1

Evidence Recovery Systems: Sexual Assault



Speculum Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- ISO17025 Compliant
- Clean Room Manufactured

The medical grade plastic ensures exceptional strength and provides safety and reliability for the patient and user. The smooth edge ensure comfort for the patient.

The product includes a swab for sampling the speculum after use, replacing the need to return the whole speculum for forensic analysis.

Recommended for use with Gelcat Lubricating Gel available from SceneSafe : K429

Code K385 Small / Pack Size 1

Code K386 Medium / Pack Size 1

Code K387 Large / Pack Size 1



Vaginal Swab Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- ISO17025 Compliant
- Clean Room Manufactured

To collect and store vaginal swab samples for subsequent forensic analysis to detect body fluids and/or DNA, lubricants or other contact traces.

Recommended to be used in conjunction with the range of SceneSafe speculum products: K385, K386, K387.

Code K384 / Pack Size 1



Proctoscope Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- ISO17025 Compliant
- Clean Room Manufactured

The medical grade plastic ensures exceptional strength and provides safety and reliability for the patient and user. The smooth edges ensure comfort for the patient. The products includes a swab for sampling the proctoscope after use, replacing the need to return the whole proctoscope for forensic analysis.

Recommended for use with Gelcat lubricating Gel available from SceneSafe: K429

Code K389 Small / Pack Size 1

Code K390 Medium / Pack Size 1



Anal Swab Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- ISO17025 Compliant
- Clean Room Manufactured

To collect and store anal swab sample for subsequent forensic analysis to detect body fluids and/or DNA, Lubricants or other contact traces.

Recommended for use with a SceneSafe Proctoscope product: K389, K390.

Code K388 / Pack Size 1

Evidence Recovery Systems: Sexual Assault



Saliva Swabbing (Spit) Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- ISO17025 Compliant
- Clean Room Manufactured

This standard kit is specifically designed for the collection of body fluids, in particular saliva from spitting incidents. Designed to deter antisocial behaviour, and to increase the chances of catching an offender, the sample can be taken from wet saliva that has been deposited on skin, clothing or other surfaces.




Medical Examination Forms Kit



- Meets FFLM Guidelines
- ISO17025 Compliant
- Clean Room Manufactured

Contains the necessary paperwork for use during the medical examination of a complainant or suspect where forensic samples are collected. Each pack contains four separate forms in triplicate.



Body Outline Forms Kit



- Meets FFLM Guidelines
- ISO17025 Compliant
- Clean Room Manufactured

Body outline diagrams, for the examiner to record observations with regard to injuries or bruising noted during a medical examination of a complainant or suspect. Contains fourteen different diagrams.



Emergency Information Kit



- Meets FFLM Guidelines
- ISO17025 Compliant
- Clean Room Manufactured

Victim information pack, which includes:

- ACPO Booklet: 'Help For Those Who Have Been Sexually Assaulted'.
- AFP 'What Happens Now?' booklet.
- BPAS 'Emergency Contraceptive Pills' leaflet.

Code K661 - Standard / Pack Size 1

Code K403 / Pack Size 1

Code K404 / Pack Size 1

Code K405 / Pack Size 1

NOT GUILTY

What is the **true cost** of contaminated evidence?

At SceneSafe we understand that the true cost of contaminated evidence is not just measured in the amount of money or man hours spent trying to conclude an investigation successfully. We know that there can be no margin of error in the identification, collection, storage and transportation of vital evidence from the scene of a crime to the laboratory, which is why we work with experts to develop a range of world-class products that meet the highest quality standard

For more information visit www.scenesafe.co.uk

Evidence Recovery Systems: Trace DNA

DNA analysis is frequently used to acquire information from biological material to aid enquiries associated with criminal offences, disaster victim identification and missing person investigations.

With the recent introduction of DNA-17, analysis techniques for DNA profiling has advanced. By using DNA-17 technology there is greater increase in potential matches from trace evidence found at a scene or on a suspect person.

Due to increased sensitivity it is paramount that stringent anti-contamination practices are followed at all times. SceneSafe follow current guidelines set in place through ISO/DIS 18385 to ensure that every step is taken to maintain evidential integrity.



Evidence Recovery Systems: Trace DNA

With the advances in DNA Profiling techniques (DNA-17) becoming increasingly more sensitive allowing forensic practitioners to obtain DNA profiles from trace DNA. Trace Recovery kits have been designed to quickly capture DNA using a quality controlled Forensic DNA Grade product. Many of our Trace DNA solutions have been developed with leading Police/Forensic units around the world.



Enhanced DNA Recovery Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- ISO17025 Compliant
- Clean Room Manufactured

For the recovery of contact DNA material. This kit contains a nylon flocked swab and special reagent that improves the recovery of DNA material from surfaces.

Code K542 / Pack Size 1



Flocked Swab DNA Recovery Kit



- Forensic DNA grade
- DNA Denatured
- Meets FFLM Guidelines
- ISO17025 Compliant
- Clean Room Manufactured

For the recovery of contact DNA material. This kit contains flocked swab, sterile water, collection tube and tamper-evident packaging.

Code K546 / Pack Size 1



Flocked Swabs



- Forensic DNA grade
- DNA Denatured
- ISO17025 Compliant

An innovative nylon flocked swab with a unique surface texture which improves the recovery of saliva and cellular material particularly from hard (non-absorbent) surfaces and therefore greatly enhance the success rates of obtaining a full profile. Items suitable for sampling would include bottles, cans, tools and other household items where the perpetrators DNA may be located.

Code K547 / Pack Size 100

Evidence Recovery Systems: Trace DNA




Scenesafe FAST™ Kit



- Forensic DNA grade
- DNA Denatured
- ISO17025 Compliant
- Clean Room Manufactured

An excellent recovery mechanism for collecting trace DNA from large surface areas, both absorbent and non-absorbent. It is ideal to be used on household items, and on items too large to submit to the laboratory. It can also be used in conjunction with swabs to maximise evidence recovery. Each kit contains five mini-lifting tapes, five DNA-free sample collection tubes, five tamper-evident return bags, instruction sheet. It can also be used in conjunction with swabs to maximise evidence recovery.

Code K543 / Pack Size 1



Scenesafe FAST™ Pack



- Forensic DNA grade
- DNA Denatured
- ISO17025 Compliant
- Clean Room Manufactured

A pack of 24 Mini-lifting tapes. Collection tubes are available to purchase separately (K480 30ml universal containers). Scenesafe FAST™ has the advantage that it does not need to be kept in frozen storage.

Code K545 / Pack Size 1

Code K544 / Pack Size 1



Contact DNA Recovery Kit



- Forensic DNA grade
- DNA Denatured
- ISO17025 Compliant
- Clean Room Manufactured

For the recovery of trace DNA material from fingerprint residues. The product has been designed to allow the scene of crime examiner to take full advantage of all available trace DNA after the fingerprint has been removed by way of aluminium powder and fingerprint lifting tape.

Code K555 / Pack Size 1

Evidence Recovery Systems: Toxicology



SceneSafe provide a comprehensive range of products to collect and submit samples for alcohol, drug and volatiles testing in relation to road traffic accidents. These kits are used where there is suspected drug or alcohol abuse and also in cases of suspicious death.

All SceneSafe toxicology products are quality control tested for the presence of alcohols, drugs of abuse, drug metabolites, medical and drug facilitated sexual assault drugs to ensure that they are free of any possible contaminants.



RTA Blood Alcohol/ Drug Kit



- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

Quality assured kits designed to enable the safe collection and transfer of a blood sample from an individual suspected of driving under the influence of alcohol and/or drugs into dual preservative coated vials (specially formulated to ensure successful analysis) providing an evidential police sample and a sample for the suspect to have privately analysed if so wished.

Code K125 - England & Wales / Pack Size 1



RTA Urine Alcohol/ Drug Kit



- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

Quality assured kits designed to enable the collection and transfer of a urine sample from an individual suspected of driving under the influence of alcohol and/or drugs into dual vials with preservative tablets (specially formulated to ensure successful analysis) providing an evidential police sample and a sample for the suspect to have privately analysed if so wished.

Code K135 - England & Wales / Pack Size 1

Evidence Recovery Systems: Toxicology



Mini Toxicology Kit



- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

A compact kit for the collection, preservation, storage and transport of toxicological body fluid samples during a post mortem examination. Preserved samples are collected in containers with optimised anti-coagulant and/or preservative to ensure a successful analysis.

Code K112 / Pack Size 1



Scenesafe™ Drugs In Hair Kit



- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

Comprehensive chain of custody hair collection kit which contains items required to collect and dispatch a hair sample to the laboratory for drugs screening.

- No special storage conditions required
- Samples collected can be screened for general drug usage or segmentation analysis.

Code K127 / Pack Size 1



Autopsy Specimen Kit



- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

A comprehensive kit for the collection, preservation, storage and transport of toxicological samples of body fluids, organs, muscle, stomach contents etc. as required during a post mortem examination. Preserved samples are collected in containers with optimised anti-coagulant and/or preservative to ensure a successful analysis.

Code K110 / Pack Size 1

Evidence Collection: Firearms

Gunshot Residue (GSR)

Primer residues can be retrieved by repeatedly sampling with adhesive tapes, or stubs, on to skin, clothing or other surfaces or by sampling suitable non-fabric surfaces with swabs. Propellant additives are recovered by either, vacuuming clothing and other surfaces, or by swabbing non-fabric surfaces. The resulting samples can also be analysed for primer.

Persistence of GSR

GSR is quickly lost from the hands but may persist for slightly longer on the face and hair. It can remain for up to a day on the other surfaces of clothing and indefinitely inside pockets depending upon the amount of disturbances.



GSR Vehicle/Scene Tape Sampling Kit



- Gunshot residue free
- ISO17025 Compliant
- Clean Room Manufactured

A tape based product for the collection of primer gunshot residue for sampling surfaces such as a vehicle or scene of crime (windows, doors etc.) or weapons. Sampling strategies for different scenarios are included in the instructions.

Code K142 / Pack Size 1



GSR Gun/Cartridge Sampling Kit



- Gunshot residue free
- ISO17025 Compliant
- Clean Room Manufactured

A 'stick' swab product for the collection of primer gunshot residue from the barrel of a gun or from discarded cartridge case(s). This evidence is useful in linking fired ammunition to a specific weapon.

Code K143 / Pack Size 1

Evidence Recovery Systems: Firearms



GSR Subject/Scene Swab Sampling Kit



- Gunshot residue free
- ISO17025 Compliant
- Clean Room Manufactured

A chemically treated swab product for the collection of primer and propellant gunshot residue from the skin and hair of an individual suspected of being involved in an incident(s) concerned with the discharge of a weapon(s) or for sampling other surfaces such as a vehicle or scene of crime (windows, doors etc.) or weapons.

Code K144 / Pack Size 1



GSR Subject Tape Sampling Kit



- Gunshot residue free
- ISO17025 Compliant
- Clean Room Manufactured

A tape based product for the collection of primer gunshot residue from the skin and hair of an individual suspected of being involved in an incident(s) concerned with the discharge of a weapon(s).
N.B. May also be used to sample individuals hands in cases of metal theft.

Code K145 / Pack Size 1



GSR Subject/Scene Stub Kit



- Gunshot residue free
- ISO17025 Compliant
- Clean Room Manufactured

A stub based product for the collection of primer gunshot residue from the skin and hair of an individual suspected of being involved in an incident(s) concerned with the discharge of a weapon(s) or for sampling other surfaces such as a vehicle or scene of crime (windows, doors etc.) or weapons. Sampling strategies for different scenarios are included in the instructions.

Code K147 / Pack Size 1

Evidence Recovery Systems: Drugs

SceneSafe offer a comprehensive range of drug detection products for the identification of drugs in hair, sweat, saliva, in powders, solids, liquids and on surfaces. Our product range includes those products approved for use under the EDIT (Evidential Drug Identification Testing) process.

Usage of any drug detection device should be carried out by trained operatives only and according to the instructions.



Evidence Recovery Systems: Drugs



Individual Drug Sampling Kit

- Drug Free
- ISO17025 Compliant



The individual sampling kits consist of a sealed foil pouch containing instructions, a customer form, a disposable microspatula, a sampling vial and a tamper-evident bag. The customer simply samples any suspect powder using just the tip of the microspatula, adds the sample to the barcoded vial, seals it in the tamper-evident bag and returns it to our forensic laboratory together with the completed customer form. We analyse the sample in our UKAS accredited laboratory and provide the customer with a report as the contents and classification of any drugs detected within 5 working days of receipt at the laboratory. Customers can opt to upgrade to fast track analysis (within 24 hours of receipt), and/or to an evidential statement rather than a report.

Code K900 / Pack Size 25



Drug Swab Testing Kit

- Drug Free
- ISO17025 Compliant



The swab testing kit is for trace drugs or areas you suspect of being used for drug taking even when no drugs can be seen. Included in this kit are instructions, a customer form, non-latex gloves, a vial of laboratory grade water, two swabs (one to be used as a blank and one for testing) and a tamper evident bag. With this kit, the customer puts the gloves on, opens the blank swab, dips this in the water and replaces the swab in its container. The user then opens the test swab, dips this in the water and swabs the area of interest before resealing the swab in its container. The user then sends both swabs (sealed in the tamper evident bag) and the customer form to us for analysis. We analyse and report as before, however there is no option for an evidential statement on the swab kit as trace evidence is rarely admissible in court.

Code K901 / Pack Size 25

Evidence Recovery Systems: Drugs



Drug Test Kit



This drug-testing kit can be used to screen samples for restricted substances such as Cannabis, LSD, Cocaine, Barbiturates, Opiates and Amphetamines.

- Drug Free
- ISO17025 Compliant

Code K898 / Pack Size 1



Marquis Test Kit



Kit contains 8 packs of 5 ampoules of Marquis Reagent, ampoule breakers, mini spatulas and reagent neutraliser to carry out testing of unknown substances to identify whether they may be Morphine, Amphetamine or Heroin. This product is approved for use under the EDIT process.

- Drug Free
- ISO17025 Compliant

Code K899 / Pack Size 1

Evidence Recovery Kits: Personal Exploitation Kit (PEK)

This is the original NATO codified PEK designed for use by NATO troops serving abroad. The PEK recovers biometric and physical evidence from suspects both living and deceased. The modular format of the kit makes recovery of evidence simple, quick and reduces the risk of contamination.

Features:

- Biometric DNA retrieval
- Fingerprint recovery
- Photographic identification
- Suspect information capture
- Physical evidence recovery
- Hair sampling



- Forensic DNA Grade
- DNA Denatured
- Drug free
- Alcohol free
- ISO17025 Compliant
- Clean Room Manufactured

We also offer a bespoke PEK which can be translated to any language.
Please speak to a member of our development team to discuss your requirements.

Code K1171 / Pack Size 1

Index

K105.....	Forensic Medical Examination Kit	Page 8
K106.....	Early Evidence Kit	Page 9
K106SW....	Early Evidence Kits with swabs	Page 9
K108.....	Early Evidence Kit (Two urine)	Page 9
K108SW....	Early Evidence Kit (Two urine) with swabs	Page 9
K110.....	Autopsy Specimen Kit	Page 18
K112.....	Mini Toxicology Kit	Page 18
K125.....	RTA Blood Alcohol/Drug Kit	Page 17
K127.....	Drugs in Hair Kit	Page 18
K135.....	RTA Urine Alcohol/ Drugs Kit	Page 17
K142.....	GSR Vehicle/Scene Tape Sampling Kit	Page 19
K143.....	GSR Gun/Cartridge Sampling Kit	Page 19
K144.....	GSR Subject/Scene Swab Sampling Kit	Page 20
K145.....	GSR Subject Tape Sampling Kit	Page 20
K147.....	GSR Subject/Scene Stub Kit	Page 20
K380.....	Alcohol/Drug Blood Kit	Page 10
K381.....	Alcohol/Drug Urine Kit	Page 10
K383.....	Mouth Collection Kit	Page 10
K384.....	Vaginal Swab Kit	Page 11
K385.....	Speculum Kit Small	Page 11
K386.....	Speculum Kit Medium	Page 11
K387.....	Speculum Kit Large	Page 11
K388.....	Anal Swab Kit	Page 11
K389.....	Proctoscope Kit Small	Page 11
K390.....	Proctoscope Kit Medium	Page 11

K392.....	Standard Swab Kit	Page 10
K392W.....	Standard Swab Kit with Sterile Water	Page 10
K403.....	Medical Examination Forms Kit	Page 12
K404.....	Body Outline Forms Kit	Page 12
K405.....	Emergency Information Kit	Page 12
K501.....	DNA Arrest/Detention Sampling Kit	Page 6
K505.....	PACE DNA Sampling Kit	Page 5
K511.....	DNA Intelligence-Led Screen Sampling Kit	Page 6
K515.....	Elimination DNA Sampling Kit	Page 5
K530.....	Police Elimination DNA Sampling Kit	Page 5
K532.....	Police Elimination DNA Sampling Kit (SPSA)	Page 6
K542.....	Enhanced DNA Recovery Kit	Page 15
K543.....	Fast Kit	Page 16
K544.....	Fast Pack	Page 16
K545.....	Fast Box	Page 16
K546.....	Flocked Swab DNA Recovery Kit	Page 15
K547.....	Flocked Swabs	Page 15
K555.....	Contact DNA Recovery Kit	Page 16
K661.....	Saliva Swabbing (Spit) Kit	Page 12
K900.....	Individual Drug Sampling Kit	Page 22
K901.....	Drug Swab Testing Kit	Page 22
K898.....	Drug Test Kit	Page 23
K899.....	Marquis Test Kit	Page 23
K1169.....	Personal Exploitation Kit (PEK)	Page 24

Terms & Conditions

How to buy:

Telephone: +44 (0)1621 786654

Fax: +44 (0)1621 786655

Email: sales@scenesafe.co.uk

SceneSafe have a network of distributors worldwide who can service your requirements, contact our sales team who can direct you to the nearest SceneSafe representative.

Shipping & Handling

For details of delivery charges outside of mainland UK please contact us. Shipping rates will be based upon the weight and dimensions of the package. Therefore charges will vary dependant on the order quantities or service required and final delivery location. Our SceneSafe freight team distribute worldwide by air, land or sea, including hazardous goods for transport.

Hazardous Goods for Transport

SceneSafe staff are IATA, IMDG & ADR certified to carry out a full Dangerous Goods documentation service to any destination specified within International shipping regulations. All hazardous products within the SceneSafe range will be identified at time of quotation.

Export Documentation

As a value added service we provide full shipping documentation required for any specific country including Certificate of Origin, Commercial Invoices, Packing Lists and Certificate of Conformity.

Terms & Conditions

For a copy of our full terms and conditions please go to our website where these can be found, alternatively contact our sales team who can provide a copy via email or fax.

scenesafe™

evidence recovery systems

Tel: +44 (0) 1621 786654

Fax: +44 (0) 1621 786655

Email: Sales@scenesafe.co.uk

www.scenesafe.co.uk

SceneSafe Ltd, registered in the UK (number 4028076), a part of the WA Group of Companies
Registered office and Headquarters: Midas House 8 & 9 Burnham Business Park, Springfield Road, Burnham-on-Crouch, Essex, CM0 8TE, United Kingdom

Authorised Distributor

